

# Easy Dental 2008 Component Update 3

Component Update 3 provides enhancements and fixes for errors reported in the latest Easy Dental release. In Component Update 3, the North Carolina NCMED3F form has been added, and errors related to printing Route Slips and surface painting in the Chart have been fixed.

For detailed information about the Easy Dental 2008 Component Update 3, please visit the Support knowledgebase at [www.easydental.com/support/knowledgebase/logon.asp](http://www.easydental.com/support/knowledgebase/logon.asp).

## Downloading and Installing Component Update 3

The component update must be installed on each workstation that runs Easy Dental 2008.

*Note: Before you install, you should have a current backup of your Easy Dental database, and Easy Dental should be closed on all workstations.*

1. Open your Internet browser and go to **www.easydental.com**.
2. Click the **Support** tab to open the Support page.
3. Click **Product Updates** in the left navigation panel to open the Product Updates page.
4. Click the **EZD2008CU3.exe** link in the Easy Dental 2008 CU 3 section of the page. The File Download dialog box appears.
5. Click **Save** to save the file to your computer. The Save As dialog box appears.
6. Expand the **Save in** drop-down list by clicking the arrow. Select **Desktop** from the list.
7. Click the **Save** button to save the file to your computer's desktop. When the file has been saved, the Download Complete dialog box appears.
8. If you want to run the update now, click the **Run** button. If a security warning appears, click **Run** to continue with the installation. Click **Yes** to install Component Update 3. Click **Finish** when the install is complete.
9. If you want to run the update later, click the **Close** button and run the update at a later time by double clicking the **EZD2008CU3** icon on your Desktop. Click the **Run** button. If a security warning appears, click **Run** to continue with the installation. Click **Yes** to install the Component Update. Click **Finish** when the install is complete.
10. Install the Component Update on every computer in your office using one of the following three methods. The first method requires an Internet connection at each computer in your office. The second option does not require an internet connection at each computer, but it requires a CD burner or external drive. The third option is more technical and should not be attempted without the help of a CIE or computer technician.

**Option 1:** Repeat steps 1 – 8 on each individual computer.

**Option 2:** Burn the Component Update file to a CD using a CD burner or copy the file to an external storage device, such as a jump drive. Take the CD or jump drive to each computer, save the EZD2008CU3 file to the Desktop, and complete steps 6 – 8 to install the Component Update.

**Option 3:** Copy the file into the EzDental folder on the server (ie. C:\Program Files\EzDental\DATA\Installs\2008). Then, from each workstation, browse to the server's EzDental Installs folder and double click EZD2008CU3.exe to install the update.

*Note: For help with the download and installation of Easy Dental Component Update 3, call Easy Dental support at (800) 824-6375, Option 2.*

## Fixes in Easy Dental 2008 Component Update 3

The North Carolina NCMED3F claim form has been added to Component Update 3. Component Update 3 also includes fixes for the following errors in Easy Dental 2008.

- Schedule:**
- Long Appointment ID numbers cause the Route Slip not to print when the appointment includes treatment-planned items.
- Chart:**
- When posting a two-surface restoration in the Chart, the restoration is painted as two one-surface restorations instead of one two-surface restoration.